SPINCOTECH SYSTEMS LLP

'Spinco Towers', 83-84, Industrial Estate, Perungudi, Chennai - 600 096. Tel: 044-4567 8947/53 Email: sales@spincotech.com www.spincotech.com LLPIN: AAI-0067 GST No.:33ADIF50121K2ZB

Ref: SLP/22/22-23 Nov 22, 2022

Parrish Pharmaceuticals

Sidcul, Haridwar

Kind Attn: Dr.Ravikant

Shimadzu P-Series High Performance All Round HPLC System with UV Detector



Dear Sir,


We Spinco represent Shimadzu Corporation, Japan exclusively in India for their LC, LCMS Systems and have installed more than 25,000 Sets of Shimadzu HPLC's, 800+ client server and 500+ LCMS/MS across the country.

Shimadzu P-Series High Performance All Round HPLC System is widely used in pharmaceutical, chemical, food and environmental analysis.

Shimadzu P-Series HPLC System is now further enhanced with additional capabilities of highest sensitive temperature controlled detectors and improved plungers and improved check-valves to deliver pulse free solvent delivery using auto pulsation correction mechanism.

It has by default heating and cooling features in its auto sampler as well as in its column oven.

Shimadzu detectors are temperature controlled to avoid any chance of base line drift due to ambient temperature fluctuation.



SPINCOTECH SYSTEMS LLP

'Spinco Towers', 83-84, Industrial Estate, Perungudi, Chennai - 600 096. Tel: 044-4567 8947/53 Email: sales@spincotech.com www.spincotech.com LLPIN: AAI-0067 GST No.:33ADIFS0121K2ZB

Shimadzu is the only company to provide a 5-Line degasser in which an extra 5th Line is provided specially for the rinsing solution which guarantee lowest RSD and lowest carryover.

ASP (Auto System Preparation) for complete Automation of system. It allows auto-startup of system, auto-purge, auto baseline stability, auto-system suitability check, auto-column equilibration and auto sample processing and auto-shut down. It allows users to use the system without their actual presence enabling unattended user operations to increase productivity of the lab.



Uniqueness of P-Series Configuration:

Pump:

- 1. The LC-20AD offers the fastest solvent delivery performance in the world. With an automatic pulsation correction mechanism and high speed micro plunger driving, it achieves pulse-free and accurate solvent delivery.
- 2. Flow Rate: 0.0001 to 10ml/min which is most precise and lowest available flow rate option for increase of column life.
- 3. Lowest Flow rate RSD of 0.06% and lowest flow rate gradient precision RSD:0.1%
- 4. Unique Parallel type double plungers design for trouble free performance.

Degasser:

- 1. A very unique 5-Line degasser which is offered only by Shimadzu.
- 2. Shimadzu provides 4+1 line degasser. 1 line extra is given to degass rinsing solution as well which helps in attaining lowest RSD and carryover and precise results.

Autosampler:

- 1. A high capacity autosampler with 175 vials capacity for 1ml vials or 105 (1.5 mL vials) etc. + 10 vials extra in any case.
- 2. Most precise injection volume accuracy with 0.3% RSD Max and injection volume accuracy of 1% max.
- 3. Lowest carryover of 0.0025% max
- 4. Both cooling and heating features used together with dehumidifying function, 4 to 40°C
- 5. Injection volume range of 0.1μl to 100μl

Column Oven:

- **1.** Both Heating as well as cooling feature in column oven with temp. Range **4-80°C**. (Room temperature -12 to 80°C).
- **2.** Slim type column oven for less space occupancy.

UV Detector:

- 1. Detector with temperature controlled mechanism (5°C above Room Temperature to 50°C of temperature control range).
- 2. Lowest noise and drift compared to any other vendor detector

'Spinco Towers', 81-84, Industrial Estate, Perungudi, Chennai - 600 096. Tel: 044-4567 8947/53 Email: sales@spincotech.com www.spincotech.com

LLPIN: AAJ-0067 GST No.:33ADIFS0121K2ZB

Ref: SLP/22/22-23 Nov 22, 2022

Parrish Pharmaceuticals Plot No -50, Haridwar

Kind Attn: Dr.Ravikant

Sub: Offer for Shimadzu P-Series HPLC system with UV Detector (1 Set)

	1	1						
Your Ref:	SLP/22/22-23		Delivery		2-3 weeks from the date of receipt of your advance payment			
Price Inco terms	FOR Destination Site		Installation & training		Will be done by trained Spinco Engineers			
Payment terms	100% Advance by RTGS Against Proforma Invoice		Warranty		12 months from the date of installation (or) 13 months from the date of invoice whichever comes earlier			
Offer Validity	10 Days		HSN	Code	9027-2000			
Shimadzu P- Series HPLC System With UV Detector Consisting Of: -								
Sr. No.	P/N	Description		QTY.	Total Price INR			
1	228-65012-58	 Pumping Method: Parallel Type Double Plunger High speed double micro-plunger parallel pump with automatic pulsation correction mechanism achieving pulse- free solvent Delivery. Flow rate range: 0.0001ml to 10.0 ml/min each Pump Built in self-diagnostic, Auto Purge function Flow precision is less than ±0.06% RSD 	1	1				
2	228-45567-91	Washing Kit, 20AD		1				
3	228-45093-93	Gradient Mixer, SUS 20A		1				
4	228-65715-41	Reservoir Tray, 20A		1				
5	228-45040-58	Quaternary Gradient Unit		1				
6	228-65555-58	DGU-20A5R BLK Degassing Unit Degassing Unit: 5 Lines (Mobile Phase 4 + Rinse Solution 1)		1				



LLPIN: AAJ-0067 GST No.:33ADIFS0121K2ZB

7	228-65126-41	Auto Sampler, SIL-20ACHT World's fastest and cleanest sampler 21 Sec injection time for 10µl sample Injection volume range-0.1ul to 100ul Injection volume precision- RSD 0.3% Near zero carryover of 0.0025% Peltier cooler programmable from 4-40C Built in dehumidifying function	1			
8	228-65205-58	CTO-10ASVP BLK Column Oven Block Heating & Electronic Cooling Method. Temperature Control Range: Room temperature −12 to 85°C.	1			
9	228-65308-58	 SPD-20A BLK UV Detector Wavelength Range: 190-700 nm Spectral bandwidth: 8nm Wavelength Accuracy: ≤±1nm Simultaneously monitoring of Dual Wavelength Linearity: 2.5 AU Noise Level: ≤±5 x 10-6 AU Drift: ≤100x10-6 AU/h Flow Cell: 12ul (10mm, TC), 12 MPa 	1			
10	223-19199-51	LabSolutions LC Connect CS Ready	1			
11	226-88585-31	Solvent Bottles with Cap 5/PKT	1			
12	228-45011-58	CBM-20Alite System Controller	1			
13	228-59156	Power Cord	4			
14	NPN	1.5 ml Vial Caps and Septa (100 Pcs/ Pkt)	1			
	TOTAL PI	INR	15,00,000			
Add IGST @18% Extra						

'Spinco Towers', 81-84, Industrial Estate, Perungudi, Chennai - 600 096. Tel: 044-4567 8947/53 Email: sales@spincotech.com www.spincotech.com

LLPIN: AAJ-0067 GST No.:33ADIFS0121K2ZB

Terms and Conditions:

Prices FOR Destination Price

2 **IGST** EXTRA @ 18% 3 Validity 30 Days

Ordering Information

Purchase Order (PO) to be released in favor of Spincotech Systems LLP Spinco Towers, 81-84, Industrial Estate, Perungudi, Chennai - 600 096.

Tamil Nadu

GST No:33ADIFS0121K2ZB

5 **Delivery** 2-3 weeks from the date of receipt of your advance payment

6 **Payment** 100% Advance through RTGS Bankers Details:

> Name of the Bank **HDFC Bank**

Name of the Account Spincotech Systems LLP

Address Old No.319B/13, New No.46B/13, New

No.46/B13 South Boag Road, T.Nagar,

Chennai - 600 017 50200024512275

Current A/C No. Account Type **Current Account** MiCR No. 600240084

IFSC Code HDFC 0002094

IMPORTANT CAUTION: Our HDFC business account 50200024512275 is a PERMANENT ACCOUNT. We DO NOT CHANGE our Bank Account used for payment remittance. Hence, if any such phishing mail reaches you seeking Bank Account change, IMMEDIATELY CALL & speak to our Order Processing Team for guidance. DO NOT create any new account/process any payments with such change of Bank Account 33ADIFS0121K2ZB

GSTN No

Order Subject to issuance of Order Acknowledgement (OA) by Spincotech Systems **Acceptance** LLP. In case of PO cancellation post OA, a 10% cancellation charge is

applicable.

Installation & Will be done by our trained and certified personnel on-site at no extra cost. Training: Our on-site visit will commence upon written confirmation of availability of

all installation prerequisites. 10 Warranty

(i) 12 months from date of installation or 13 months from the date of

Invoice, whichever is earlier.

(ii) Shimadzu/Spinco will not be liable for system performance and any

malfunctioning caused by usage of unauthorized third-party

consumables/repair parts. In such cases, the system warranty will become

null and void.

Yours faithfully,

for SPINCOTECH SYSTEMS LLP

Ayush Singh

Sr.Sales Executive

avush.singh@spincotech.com

Mobile: 7080432553